


Disposable Biopsy Forceps

Indications

It is used for endoscopic sampling of living tissue and removal of foreign bodies in the digestive tract and respiratory tract.



Bite Block



Operating Sheet

Features

- ※ The outer surface is coated with a special material super-smooth coating, which can be inserted smoothly into the forceps channel, effectively reducing the wear of the instrument on the endoscope forceps channel.
- ※ Internal flexible connection can pass through any curved endoscope channel
- ※ Unique four-bar plier head design, the jaws are firmly engaged and not easy to slip off
- ※ Ergonomic handle design provides excellent operating feel
- ※ Sterile kit packaging, easy to use



Plain cup forceps



Plain cup forceps with needle



Serrated forceps



Serrated forceps with needle



Specifications

Code	Disposable Biopsy Forceps				Bite Blocks (mm) h × d × l	Operating Sheet (cm)
	Working Length (mm)	Cup (Jaw) O.D. (mm)	Cup Opening Size (mm)	Cup Length (mm)		
BF18-B18/120P	1200 ± 50	1.8 ± 0.5	5	2	26 × 22 × 25 ± 2	35 × 45 ± 2
BF18-G23/160P	1600 ± 50	2.3 ± 0.5	6.5	3.5		
BF18-G23/180P	1800 ± 50	2.3 ± 0.5	6.5	3.5		
BF18-E23/230P	2300 ± 50	2.3 ± 0.5	6.5	3.5	—	60 × 50
BF18-C23/50P	500 ± 50	2.3 ± 0.5	6.5	3.5		40 × 50 ± 5