

Disposable Polypectomy Snare


Indications

This product is used for the excision of gastrointestinal polyps under the endoscope in clinical practice.



Features

- ※ Rotation function improves polyp retrieval efficiency
- ※ A variety of different shapes and different specifications are suitable for different clinical needs
- ※ Using nickel-titanium alloy wire, it is suitable for the treatment of complex location lesions, and it is not deformed after being bent for many times.
- ※ Coil hardness is moderate, recovery is good, easy to open and close



Specifications

Loop Shape	Code	Sheath O.D.	Loop WidthA	Working Length(mm)
		±0.5mm	± 3mm	± 100mm
Ellipse	ES18-E1810/2400	Φ 1.8	10	2400
	ES18-E1815/2400	Φ 1.8	15	2400
	ES18-E1810/1800	Φ 1.8	10	1800
	ES18-E1815/1800	Φ 1.8	15	1800
	ES18-E2410/2400	Φ 2.4	10	2400
	ES18-E2415/2400	Φ 2.4	15	2400
	ES18-E2410/1800	Φ 2.4	10	1800
	ES18-E2415/1800	Φ 2.4	15	1800
	ES18-E2425/1600	Φ 2.4	25	1600
	ES18-E2425/2400	Φ 2.4	25	2400
	ES18-E2432/2400	Φ 2.4	32	2400
	ES18-E2425/1800	Φ 2.4	25	1800
ES18-E2432/1800	Φ 2.4	32	1800	
Hexagonal	ES18-H2410/2400	Φ 2.4	10	2400
	ES18-H2415/2400	Φ 2.4	15	2400
	ES18-H2422/2400	Φ 2.4	22	2400
	ES18-H2432/2400	Φ 2.4	32	2400
	ES18-H2410/1800	Φ 2.4	10	1800
	ES18-H2415/1800	Φ 2.4	15	1800
	ES18-H2422/1800	Φ 2.4	22	1800
	ES18-H2432/1800	Φ 2.4	32	1800
Crescent	ES18-C2415/2400	Φ 2.4	15	2400
	ES18-C2425/2400	Φ 2.4	25	2400
	ES18-C2432/2400	Φ 2.4	32	2400
	ES18-C2415/1800	Φ 2.4	15	1800
	ES18-C2425/1800	Φ 2.4	25	1800
	ES18-C2432/1800	Φ 2.4	32	1800
Round	ES18-R2432/2400	Φ 2.4	32	2400
	ES18-R2432/1800	Φ 2.4	32	1800